

Adapter Selection Guide

SCHOTT stereo microscopy illumination



Distribution in the UK & Ireland



Lambda Photometrics Limited
Lambda House Batford Mill
Harpenden Herts AL5 5BZ
United Kingdom
E: info@lambdaphoto.co.uk
W: www.lambdaphoto.co.uk
T: +44 (0)1582 764334
F: +44 (0)1582 712084

SCHOTT
glass made of ideas

Find the right SCHOTT adapter for your microscope



Optimal illumination of objects under examination is one of the most decisive factors for good results in stereo microscopy. Ideal light is necessary for detection of hidden details, differentiation of specific areas and emphasis of special characteristics.

In order to make your selection of suitable lighting a little easier, we have put together this assortment of auxiliary templates. Print it out in actual size and find the right adapter.

Please fold here to use as a ruler (in cm)



Adapters for EasyLED and VisiLED Transmitted Light Stages

Hold your white stage plate in the middle of the copy to measure its size. Select the product number for the required adapter accordingly.

Select actual size in
your printer settings
and print it out!

Ø 180 mm

Ø 120 mm

Ø 100 mm

Ø 90 mm



Ø 180 mm | P/N 157.620

Ø 120 mm | P/N 157.610

Ø 100 mm | P/N 157.600


Ø 90 mm | P/N 157.630

SCHOTT
glass made of ideas

Adapters for Ringlights


Hold the different diameters to your lens or use the ruler to measure its size.
Select the product number for the required adapter accordingly.

Select actual size in
your printer settings
and print it out!




Ø 66 mm

Ø 66 mm | No adapter needed




Ø 60 mm

Ø 60 mm | P/N 157.452



Ø 58 mm

Ø 58 mm | P/N 157.450



Ø 50 mm

Ø 50 mm | P/N 157.448


Please fold here to use as a ruler (in cm)



SCHOTT
glass made of ideas

Adapters for KL 300 LED and EasyLED Spotlight Plus

Hold the different diameters to your column or use the ruler to measure its size.
Select the product number for the required adapter accordingly.




Ø 35 mm | P/N 120.235



Ø 32 mm | P/N 120.232

Select actual size in
your printer settings
and print it out!



Ø 29 mm | P/N 120.229



Ø 25 mm | P/N 120.225



Ø 20 mm | P/N 120.220

Please fold here to use as a ruler (in cm)



Distribution in the UK & Ireland



**Characterisation,
Measurement &
Analysis**

Lambda Photometrics Limited
Lambda House Batford Mill
Harpenden Herts AL5 5BZ
United Kingdom
E: info@lambdaphoto.co.uk
W: www.lambdaphoto.co.uk
T: +44 (0)1582 764334
F: +44 (0)1582 712084



If your microscope does not have a round
column, please contact your microscope
manufacturer to find the right adapter.

SCHOTT
glass made of ideas